


# PHARMACY FRIDGE

## USERS INSTRUCTION MANUAL




Thank you for buying this pharmacy fridge. It has been designed and manufactured for many years of trouble-free service.

Operation is very simple, but please take a few minutes to read this booklet. It contains safely information, and tips to ensure you receive from your fridge.

### Portage and installation

- When carrying the fridge , keep the fridge vertical, the angle of inclination between sides of the fridge and horizontal plane should be over 50° , otherwise it will cause trouble of the compressor and effect normal operation of the fridge .
- After unpacking the fridge, remove the packing materials and place the fridge on a level and sturdy floor near power socket.
- Do not place your fridge near a direct heat source or direct sunlight and leave at least 6cm all round the unit for adequate ventilation (wall mounted models do not need extra space at the back because the cooling system is in the sides). Failure to do so may effect the performance of your fridge.
- Never place the fridge on a location where water can reach or with high moisture.
- Ensure your fridge is level by adjusting the feet at the bottom of your fridge.

### Exterior view



## Power supply preparation


The fridge must use a special three-hole socket with fuse and well grounded. never ground the fridge through a gas pipe, tap water pipe etc. Never share same circuit with other appliances. Before plugging in your fridge for the first time make sure the power cable is free from damage and check that your voltage corresponds to that on the serial plate.

## Plug in

Once plugged in the display will show the temperature inside the fridge and the compressor will run late on (the indicator light of runing will on at the same time ). This temperature is the actual temperature recorded by the probes inside the fridge and as the fridge cools to the set point (factory set at 4°C/39 ) so the number on the display will descend accordingly.

If power is cut off, wait for at least 5 minutes before plugging the unit in again to avoid damaging the compressor.

## Front panel display



## Adjust operating temperature

Press **K4** button and hold for more than 3 seconds, the current set point temperature will flash display , and then press **K1** button or **K2** button to change the set point, you can select the desired target temperature and then stop operating any button, three seconds later the display will return to normal.



The display shows set point temperature normally, press **K3** button once the temperature inside the fridge will be showed.

### High/Low Temperature Alarm

If the temperature in your fridge rises above 9°C/48°F for more than about 7 minutes an alarm will sound and the front panel display will flash **HI** and the current temperature inside the fridge. This information will display alternately. (Plug in the fridge for the first time, If the temperature in fridge is above 9°C/48°F alarm will appear immediately.)


The cause of the rise in temperature should be investigated immediately. Usually, it is simply because the door has been open for a long time or the fridge has been restocked, however, it is recommended that you check the fridge to make sure that the temperature is returning to normal after the alarm has been activated.

If you want to silence the alarm, simply press **K1** or **K2** button.



If the temperature in your fridge drops below 1°C/34°F for more than about 7 minutes an alarm will sound and the front panel display will flash **LO** and the current temperature will be displayed alternately. There could be a number of reasons for this, such as extremely cold items being placed in the fridge.



The cause of the drop in temperature should be investigated immediately and the fridge should be monitored to check that it is returning to the correct temperature.



## Min/Max Temperature Recording

The fridge records the maximum high and minimum low temperatures when connected to mains electricity. This information is stored until the fridge is reset.

This information can be accessed if you press K3 button and release. The display will show **HI** and then the latest 30 minutes high temperature recorded since the last reset, then **LO** and display the low temperature.




## Min/Max Temperature Reset

Press **K3** button and hold for more than 3 seconds, An alarm will sound and then release **K3** button, This means you have successfully reset the minimum/maximum temperature display and the previous high and low temperatures have been deleted.

## Door Alarm

Open the door of the fridge and then the door indicator light will be on, if the door is left open for more than 4 minutes an audible alarm will sound and the front panel will flash **do**. If you want to silence the alarm, for example, you are stocking the fridge, simply press **K3** button or **K4** button.

Close the door, the door alarm will be eliminated automatically.



## Battery

Insert 4 AA 1.5V long life alkaline batteries (not supplied) in the holder located on the right side of the unit.

When power supply is cut off batteries will power the controller.


Note: batteries will **only** power the controller not the fridge.

## Remote alarm

In some cases of alarm, for example, high temperature alarm, low temperature alarm, door alarm, etc, all alarm will output to the connector located on the back of the unit.

The maximum rating for the connection is 1A 30Vdc.

**DO NOT CONNECT TO MAINS VOLTAGE.**



## Internal light

You can turn the light on or off with the switch located at the front panel,

## Storage of articles

- The fridge is not overloaded. It is recommended the fridge is only partially filled (around 55%) and even at maximum stock level the content should not exceed 70% of the fridges capacity.
- Leave some space around the goods to let cold air inside the cabinet flow. Ensure the articles do not touch the back of the interior so as not to freeze.
- Store the articles away from the inner fan.
- Do not put warm goods into the fridge directly until they are cooled to normal temperature.
- Shelves can be adjusted according to different applications.

## **Maintenance**

**It is essential that you unplug the fridge from the mains and transfer the contents where they can be stored and monitored at the correct temperatures before cleaning.**

- Wash the inner compartment with warm water and neutral detergent. DO NOT allow the control panel, cables or plug to get wet. NEVER use corrosive detergents, wire brushes, or abrasive scourers to clean your fridge. NEVER use metal or sharp implements to remove debris
- Dry all surface thoroughly.
- To ensure trouble-free operation the condenser should be cleaned every three months where appropriate using a vacuum hose. The condenser is located at back of the cabinet. In exceptionally dusty locations the condenser should be cleaned more often.

### **WARNING:**

- Do not damage the refrigerant circuit.
- Do not allow children to play with the appliance, or to sit on it or to hang onto the door.
- For safety, the appliance must be properly earthed in accordance with specifications.
- Do not damage any parts of the appliance which carry refrigerant by piercing, perforating, crushing, twisting or scraping. If the refrigerant comes into contact with the eyes, it may cause serious eye injury
- Always remember to unplug the appliance before cleaning. Never unplug your appliance by pulling on the power cable. Always grip plug firmly and pull straight out from the socket. Always check that the plug and cable are undamaged.
- Do not use electrical appliances inside the fridge, unless they are of the type recommended by the manufacturer.
- If the supply cord is damaged, it must be replaced by a special cord or assembly available from the manufacturer or its service agent.
- This appliance is not intended for use by persons (including children) with reduced physical, sensory or mental capabilities, or lack of experience and knowledge, unless they have been given supervision or instruction concerning use of the appliance by a person responsible for their safety. Children should be supervised to ensure that they do not play with the appliance.
- All repairs must be carried out by qualified engineers. Inadequate repairs may be dangerous. Never try to repair the fridge yourself.
- Do not allow children to play the packaging, and destroy plastic bags safely.

## Troubleshooting guide

Problems	Possible causes	Ways to solve problem
The unit does not operate	Bad connection of plug or fuse burn out	Unplug the power cord and re-connect the plug to power supply or replace with a new fuse
Abnormal operation of the compressor or has buzz sound	The power voltage is out of rated range	Disconnect the power supply immediately and reconnect with it after normal voltage. It is necessary to have a voltage stabilizer in case of poor power supply.
Compressor operates for a long time and not frost on the surface of the evaporator	Refrigeration system is at fault (leakage or blocked)	Call for service
There is frost or ice on the walls of the inner cabinet and internal temperature is too low, the compressor never stops running.	Thermostat does not work	Call for service
	The temperature setting of the thermostat is too low.	Adjust to a higher temperature
The internal temperature is too high and the compressor never stops operating.	Bad heat dissipation and ventilation of condenser	Improve ventilation
	Too much warm goods was put in at one time	Remove some goods so air can circulate
	Door is being opened too frequently during initial coll down	Permit the unit to cool down adequately, prevent product access during this phase
Too noisy	The unit is not level	Adjust the adjustable legs
	The fastener of the unit is loose	Fasten the loose fastener
	Pipe near the compressor touch	Carefully separate the touched pipes
The side of fridge is hot	The condenser in the side wall gives out heat as part of normal operation	nothing to worry about.
Sometimes a light sound of water flowing will be heard	Refrigerant flowing inside the pipe.	nothing to worry about.
There may be condensation on the glass door	a high ambient temperature or humid conditions,	Dried with a cloth
The display shows "OC" and compressor does not run	The temperature sensor is open circuit	Call for service
The display shows "SC" and compressor does not run	The temperature sensor is short circuit	Call for service