



**NULINE**  
REFRIGERATION

A SIGN OF QUALITY

## HF Spark Safe



Nuline Refrigeration's HF **Spark Safe** Medical & Laboratory is the safest way to store your product. Static refrigerated shelves help keep even temperature throughout. Standard with door lock for added peace of mind.

### Features & Benefits

- Audible alarm
- Static Refrigerated shelves
- Door lock
- Digital Temperature indicator

### Optional

- Chart recorder
- Data logger
- Auto dialler
- Voltage free connection



Model	Dimensions	Shelves	Capacity	Door
HFI 200	600x620x850	135	3	Solid
HFI 400	600x620x1900	350	5	Solid
HFI 600	770x750x1900	570	5	Solid

\*Nuline Refrigeration reserves the right to alter specifications without prior notice.