

The ENLAKE logo consists of the word "ENLAKE" in a bold, black, sans-serif font, centered within a white rectangular box. This box is set against a larger, orange oval background.

# Medical Refrigeration Solutions

Catalogue

5 Deakin Crt  
Wodonga VIC 3690

A smaller version of the ENLAKE logo, featuring the word "ENLAKE" in black text on a white background, enclosed in an orange oval.

P: (02) 6059 1438

F: (02) 6059 8476

E: [admin@enlake.com.au](mailto:admin@enlake.com.au)

W: [www.enlake.com.au](http://www.enlake.com.au)





## Medical Refrigeration

For vaccine and pharmaceutical storage in hospitals, pharmacies and medical practices

- under counter
- upright
- multi-door
- breast milk
- clinical trials
- custom made
- pharmacy freezers

## Laboratory Refrigerators and Freezers

For storage of temperature sensitive laboratory products and flammable liquids

- fan forced
- static
- intrinsically safe
- custom made
- stainless steel

ENLAKE™

ENLAKE are the medical and laboratory refrigeration specialists. If you don't see the product that you are looking for, please contact us as we can option and customise fridges and freezers to meet just about any need and this includes complete custom made fridges and freezers.

5 Deakin Crt  
Wodonga VIC 3690

P: (02) 6059 1438  
F: (02) 6059 8476  
E: [admin@enlake.com.au](mailto:admin@enlake.com.au)  
W: [www.enlake.com.au](http://www.enlake.com.au)



## Under Counter Pharmacy/Medical/Vaccine Refrigerators

These fridges can be conveniently placed under a counter or on top of a bench.

- ICS Pacific G135L
- Nuline HR200 & HR200G
- Liebherr LKUexb 1610 & LKUexb 1612
- Quirks HB80



5 Deakin Crt  
Wodonga VIC 3690

P: (02) 6059 1438  
F: (02) 6059 8476  
E: [admin@enlake.com.au](mailto:admin@enlake.com.au)  
W: [www.enlake.com.au](http://www.enlake.com.au)



## ICS PACIFIC G135L MEDICAL Display Refrigerator

### Features

- Full length internal rear wall grill
- Raised bottom shelf
- Superior air flow throughout the cabinet
- Air circulation fan in compressor compartment
- Internal illumination
- High/low temperature display with reset button
- Electronic temperature control and display
- Audible alarm system
- External alarm facility
- Door lock
- Self-closing door
- Maintenance free condenser
- Ventilated cooling system
- Automatic defrost
- Four adjustable shelves plus base
- Three years parts and labour warranty\* (\*Conditions apply)

### Dimensions

Width- 56cm, Depth - 62cm, Height - 85cm

### Capacity

135 litres

### Temperature Range

+2°C to +8°C



# HR Pharmacy



## NULINE REFRIGERATION HR200 and HR200G Vaccine Fridges

Nuline Refrigeration's HR Series Pharmaceutical Refrigerators have been tested and approved by The Pharmacy Guild of Australia to meet every requirement and to ensure peace of mind. All refrigerators are fitted with a controller that enables you to check your temperature min/max at any stage, simply by pushing a button.

### Features

- QCCP approved
- Digital Temperature indicator
- Built in alarm
- Door lock
- 4 shelves

### Optional

- Glass door – HR200G / Solid door – HR200
- Data logger
- Voltage Free Contacts – for connection to a BMS

### Dimensions

Width- 60mm, Depth – 620mm, Height – 850mm

### Capacity

135 litres

### Temperature Range

+2°C to +8°C



HR200



HR200G

**ENLAKE**

5 Deakin Crt  
Wodonga VIC 3690

P: (02) 6059 1438  
F: (02) 6059 8476  
E: [admin@enlake.com.au](mailto:admin@enlake.com.au)  
W: [www.enlake.com.au](http://www.enlake.com.au)

# LIEBHERR

Wherever floor space is limited or under-worktop installation is planned, the Austrian made Mediline medical refrigerators with Comfort electronic controllers are the ideal solution. These meet the exacting demands made of medical appliances.

The temperature can be set from +3°C to +8°C. The forced-air cooling system in conjunction with the precision electronic controller ensures temperature consistency and uniform temperature distribution in the interior. The new



## LIEBHERR - LKUv 1610 LKUv 1612

### Features

**Precision Comfort electronic controller with:**

**Digital temperature display**

**Integrated alarm systems – visual & audible**

**Integrated data memory with min/max temperatures**

**External alarm facilities**

**Maximum temperature stability and temperature consistency**

**Internal illumination**

**Self-closing doors with an integrated door lock**

**Reversible door**

**Automatic defrost**

**Three fully adjustable shelves + base**

**24 months warranty**

### Dimensions

**Width- 600mm, Depth – 615mm, Height – 820mm**

### Capacity

**141 litres**

### Temperature Range

**+3°C to +8°C**



5 Deakin Crt  
Wodonga VIC 3690  
P: (02) 6059 1438  
F: (02) 6059 8476  
E: [admin@enlake.com.au](mailto:admin@enlake.com.au)  
W: [www.enlake.com.au](http://www.enlake.com.au)

# QUIRKS *Medisafe*



## QCPP Vaccine Fridges G2 Model HB80

### Description

New and improved Medisafe G2 is now available. Fridges specifically designed for vaccine storage. Temperature control system preconfigured to hold temperature from 2°C to 8°C (approved QCPP specification). Easily download temperature history to any computer via the data logger.

### Features

- QCPP accredited vaccine fridge
- Electronic temperature control and monitoring system
- Audible auto resetting alarms
- Min/Max Temperature Recorder
- Easy to use Computer Data Logger
- Auto Dialer (option)
- Australian design
- 24 months warranty
- Dedicated national service and help line  
1800 686 672

### Dimensions

Width- 50mm, Depth – 560mm, Height – 830mm

### Capacity

80 litres

### Temperature Range

+2°C to +8°C



ENLAKE

5 Deakin Crt  
Wodonga VIC 3690

P: (02) 6059 1438  
F: (02) 6059 8476  
E: [admin@enlake.com.au](mailto:admin@enlake.com.au)  
W: [www.enlake.com.au](http://www.enlake.com.au)



## Upright Pharmacy/Medical/Vaccine Refrigerators

These are the most commonly used vaccine fridges, all of which are approved by the Pharmacy Guild of Australia.

- ICS Pacific G315L
- Nuline HR400 & HR400G
- Nuline HR600 & HR600G
- Liebherr LKv 3910 & LKv 3912
- Quirks FKG311 & FKG371



5 Deakin Crt  
Wodonga VIC 3690

P: (02) 6059 1438  
F: (02) 6059 8476  
E: [admin@enlake.com.au](mailto:admin@enlake.com.au)  
W: [www.enlake.com.au](http://www.enlake.com.au)





## ICS PACIFIC G315L MEDICAL Display Refrigerator

### Features

- 100mm insulation = energy consumption of only 3.1kw/24hr
  - LED illumination set into the glass door
  - High/low temperature display with reset button
  - Electronic temperature control and display
  - Audible alarm system
  - External alarm facility
  - Glass or solid door able to be hinged on either side
  - Door lock
  - Efficient ventilated cooling system
  - Automatic defrost
  - 7 adjustable shelves and 3 small shelves including base
  - Maintenance free condenser
- Three years parts and labour warranty\***  
(\*Conditions apply)

### Dimensions

Width- 59cm, Depth - 64cm, Height - 187cm

### Capacity

315 litres

### Temperature Range

+2°C to +8°C

**NEW ENERGY EFFICIENT MODEL**



5 Deakin Crt  
Wodonga VIC 3690

P: (02) 6059 1438  
F: (02) 6059 8476  
E: [admin@enlake.com.au](mailto:admin@enlake.com.au)  
W: [www.enlake.com.au](http://www.enlake.com.au)



# HR Pharmacy



HR400

**ENLAKE**

5 Deakin Crt  
Wodonga VIC 3690

P: (02) 6059 1438  
F: (02) 6059 8476  
E: [admin@enlake.com.au](mailto:admin@enlake.com.au)  
W: [www.enlake.com.au](http://www.enlake.com.au)

## **NULINE REFRIGERATION HR400 and HR400G Vaccine Fridges**

Nuline Refrigeration's HR Series Pharmaceutical Refrigerators have been tested and approved by The Pharmacy Guild of Australia to meet every requirement and to ensure peace of mind. All refrigerators are fitted with a controller that enables you to check your temperature min/max at any stage, simply by pushing a button.

### **Features**

- QCCP approved
- Digital Temperature indicator
- Built in alarm
- Door lock
- 6 shelves

### **Optional**

- Glass door – HR400G / Solid door – HR400
- Data logger
- Voltage Free Contacts – for connection to a BMS

### **Dimensions**

Width- 600mm, Depth – 620mm, Height – 1850mm

### **Capacity**

350 litres

### **Temperature Range**

+2°C to +8°C



HR400G

## NULINE REFRIGERATION HR600 and HR600G Vaccine Fridges

Nuline Refrigeration's HR Series Pharmaceutical Refrigerators have been tested and approved by The Pharmacy Guild of Australia to meet every requirement and to ensure peace of mind. All refrigerators are fitted with a controller that enables you to check your temperature min/max at any stage, simply by pushing a button.

### Features

- QCCP approved
- Digital Temperature indicator
- Built in alarm
- Door lock
- 6 shelves

### Optional

- Glass door – HR600G / Solid door – HR600
- Data logger
- Voltage Free Contacts – for connection to a BMS
- Castors

### Dimensions

Width- 770mm, Depth – 750mm, Height – 1900mm

### Capacity

570 litres

### Temperature Range

+2°C to +8°C



HR400G



5 Deakin Crt  
Wodonga VIC 3690  
P: (02) 6059 1438  
F: (02) 6059 8476  
E: [admin@enlake.com.au](mailto:admin@enlake.com.au)  
W: [www.enlake.com.au](http://www.enlake.com.au)

# LIEBHERR

LKUv3912



LKUv3910

## LKUv 3910 LKUv 3912

### Features

**Precision Comfort electronic controller with:**  
**Digital temperature display**  
**Integrated alarm visual & audible alarms**  
**Integrated data memory**  
**Min/max temperature recording**  
**External alarm facilities**

**Maximum temperature stability & consistency**  
**Internal illumination**  
**Self-closing doors with an integrated door lock**  
**Reversible door**  
**Automatic defrost**  
**Five fully adjustable shelves + base**  
**24 months warranty**

### Dimensions

**Width- 600mm, Depth – 615mm, Height – 1840mm**

### Capacity

**T**he Austrian made Mediline medical refrigerators with Comfort electronic controllers is the ideal vaccine fridge.

The temperature can be set from +3°C to +8°C. The forced-air cooling system in conjunction with the precision electronic controller ensures temperature consistency and uniform temperature distribution in the interior. The new

**M**aximum temperature stability and consistency. The forced-air cooling system in conjunction with the precision Comfort electronic controller ensures maximum temperature stability and high temperature consistency in the interior. The appliances are with a secondary safety thermostat to prevent temperatures below +2°C in case of a fault. Tested according to the standard with regard to maximum temperature stability and optimal temperature consistency.



5 Deakin Crt  
Wodonga VIC 3690  
P: (02) 6059 1438  
F: (02) 6059 8476  
E: [admin@enlake.com.au](mailto:admin@enlake.com.au)  
W: [www.enlake.com.au](http://www.enlake.com.au)

# QUIRKS Medisafe



## QCPP Vaccine Fridges G2 Models 311, 371

### Description

New and improved Medisafe G2 is now available. Fridges specifically designed for vaccine storage. Temperature control system preconfigured to hold temperature from 2°C to 8°C (approved QCPP specification). Easily download temperature history to any computer via the data logger.

### Features

- QCPP accredited vaccine fridge
- Electronic temperature control and monitoring system
- Audible auto resetting alarms
- Min/Max Temperature Recorder
- Easy to use Computer Data Logger
- Auto Dialer (option)
- Australian design
- 24 months warranty
- 311 & 371 comes standard with lock
- Dedicated national service and help line 1800 686 672

### Dimensions

Width- 595mm, Depth – 595mm  
Height –1630mm (311) 1930mm (371)

### Capacity

306 litres (311), 381 litres (371)

### Temperature Range

+2°C to +8°C

311



312



**ENLAKE**

5 Deakin Crt  
Wodonga VIC 3690

P: (02) 6059 1438  
F: (02) 6059 8476  
E: [admin@enlake.com.au](mailto:admin@enlake.com.au)  
W: [www.enlake.com.au](http://www.enlake.com.au)