



ENLAKE Schedule 8 Drug Refrigerator

This refrigerator is designed to meet NSW Health's new guidelines for storing Schedule 8 (S8) drugs.

Features specific to the storage of Schedule 8 Drugs

- This refrigerator has brackets that allow it to be bolted to the floor and the wall thus allowing it to be securely attached to the premises (bolts and labour are not included)
- A Schedule 8 drug fridge is to be securely locked when not in use. For this reason, the door will automatically lock as soon as it closes. This prevents the fridge from ever being inadvertently left unlocked.
- It is critical that access to Schedule 8 Drug Refrigerator is limited to only authorised people such as a specific pharmacist, nurse, midwife or medical practitioner. The normal locks found on a pharmacy fridge are not able to satisfactorily meet this requirement as most pharmacy fridge manufacturers key their entire product range alike. This means that there may be other keys on site in the hands of unauthorised people that may be able to access the fridge. In addition, such keys are not protected against unauthorised key duplication. It is for this reason that our fridge is unlocked by your choice of a keypad or a swipe card.
- A solid door keeps the fridge contents hidden from prying eyes



Lock Features

- Keep unauthorised personnel from accessing your pharmacy refrigerator
- No keys, padlocks or loose cables
- Silent and failsafe operation makes it easy to use
- Only a four digit code or swipe needed for easy access
- Uses a magnetic lock with no moving parts to wear
- The Keypad can be programmed with different access codes and this can be changed as often as you like
- Manufactured from heavy duty zinc, alloy electroplated, which makes it vandal resistant and moisture proof



Features related to the storage of medication

- Quality Care Pharmacy Program (QCPP) approved
- TGA approved
- Holds the product within 2 to 8°C
- Built in alarms
- Digital temperature indicator
- Fan forced for even temperature
- Door lock
- 12 months warranty for parts and labour
- RC5 USB Data logger

Specifications

Temperature range	2–4°C
Capacity (litres)	125
Refrigerant	R134a
Shelves	4
Temperature indicator	Digital
Unit dimensions (WxDxH)	600 x 600 x 825 (width excludes swipe/keypad)
Door	Solid

